Literature references

- 1 Schwarz et al. Journal of Clinical Periodontology. 2006;33:491-499. Healing of intrabony peri-implantitis defects following application of a nanocrystalline hydroxyapatite (Ostim) or a bovine-derived xenograft (Bio-Oss) in combination with a collagen membrane (Bio-Gide). A case series.
- 2 Schwarz et al. Journal of Clinical Periodontology. 2006;33:584-595. Influence of different treatment approaches on non-submerged and submerged healing of ligature induced peri-implantitis lesions: an experimental study in dogs.
- 3 Schwarz et al. Clinical Oral Implants Research. 2007;18:161-170. Comparison of naturally occurring and ligature-induced peri-implantitis bone defects in humans and dogs.
- 4 Schwarz et al. Journal of Clinical Periodontology. 2008;35:80-87. Two-year clinical results following treatment of peri-implantitis lesions using a nanocrystalline hydroxyapatite or a natural bone mineral in combination with a collagen membrane.
- 5 Schwarz et al. Journal of Clinical Periodontology. 2009;36:807-814. Surgical regenerative treatment of peri-implantitis lesions using a nanocrystalline hydroxyapatite or a natural bone mineral in combination with a collagen membrane: a four-year clinical follow-up report.
- 6 Schwarz et al. Journal of Clinical Periodontology. 2010;37:449-455. Impact of defect configuration on the clinical outcome following surgical regenerative therapy of periimplantitis.
- 7 Schwarz & Becker. Quintessence Publishing Co. Ltd. London, 2010. Peri-implant Infection. Etiology, Diagnosis and Treatment.

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Indication Sheet

Failure Management

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Treatment concepts for Peri-Implantitis by PD Dr. Frank Schwarz, Dr. Narja Sahm and Prof. Dr. Juergen Becker, Duesseldorf, Germany

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
2. Synopsis: Surgical Treatment of Peri-implantitis^{3,7}

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


Surgical Procedure: Implantoplasty + Regenerative Treatment (Class Ib)^{1,4,5,6}

Surgical Procedure: Implantoplasty + Regenerative Treatment (Class Ic)^{1,4,5,6}

Soft Tissue Management: Collagen Matrix (Class Ib + II)^{1,4,5}

Surgical Procedure: Implantoplasty + Regenerative Treatment (Control) (Class Ie + II)^{1,4,5}

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©Geistlich Pharma AG