## Literature

- <sup>1</sup> Buser D, Ingimarsson S, Dula K, Lussi A, Hirt HP, Belser UC. Long-term stability of osseointegrated implants in augmented bone: A 5-year prospective study in partially edentulous patients. *Int J Periodontics Restorative Dent* 2002; 22:109-117.
- <sup>2</sup> Urban IA, Jovanovic S, Lozada JL. Vertical ridge augmentation using guided bone regeneration (GBR) in three clinical scenarios prior to implant placement: A retrospective study of 35 patients 12 to 72 months after loading. *Int J Oral Maxillofac Implants* 2009;24:502-510.
- <sup>3</sup> Zitzmann NU, Schärer P, Marinello CP. Long-term results of implants treated with guided bone regeneration: A 5-year prospective study. *Int J Oral Maxillofac Implants* 2001;16:355-366.
- <sup>4</sup> Hämmerle CHF, Jung RE, Yaman D, Lang NP. Ridge augmentation by applying bioresorbable membranes and deproteinized bovine bone mineral: A report of twelve consecutive cases. *Clin Oral Impl Res* 2008;19:19-25.
- <sup>5</sup> Urban IA, Nagursky H, Lozada JL. Horizontal ridge augmentation with a resorbable membrane and particulated autogenous bone with or without anorganic bovine bone-derived mineral: A prospective case series in 22 patients, *Int J Oral Maxillofac Implants* 2011;26(2):404-14.
- <sup>6</sup> Rothamel D, Schwarz F, Sculean A, Herten M, Scherbaum W, Becker J. Biocompatibility of various collagen membranes in cultures of human PDL fibroblasts and human osteoblast-like cells. *Clin Oral Implants Res*. 2004;15(4):443-9.
- <sup>7</sup> Schwarz F, Rothamel D, Herten M, Wüstefeld M, Sager M, Ferrari D, Becker J. Immunohistochemical characterization of guided bone regeneration at a dehiscence-type defect using different barrier membranes: An experimental study in dogs. *Clin Oral Implants Res*. 2008;19(4):402-15.

## Suppliers

Anti-inflammatory medication: 50 mg diclofenac, Cataflam®, Novartis Pharmaceuticals

Local anesthetic: Artican-hydrochloride with adrenaline 1/100,000

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
## Indication Sheet

# Horizontal Augmentation

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## Treatment concept of Dr. Istvan Urban, Loma Linda University, USA

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## 1. Indication profile

|                       |   |   |
|-----------------------|---|---|
| Region                | <input type="checkbox"/> aesthetic region                               | <input checked="" type="checkbox"/> non-aesthetic region                        |
|                       | <input type="checkbox"/> single tooth gap                               | <input type="checkbox"/> multiple tooth gap                                     |
| Bony situation        | <input checked="" type="checkbox"/> bone defect present                 | <input type="checkbox"/> no bone defect present                                 |
| Soft tissue situation | <input type="checkbox"/> recession                                      | <input type="checkbox"/> no recession   |
|                       | <input type="checkbox"/> inflamed                                       | <input type="checkbox"/> infected   |
|                       | <input type="checkbox"/> thick biotype                                  | <input type="checkbox"/> thin biotype   |
|                       | <input checked="" type="checkbox"/> primary wound closure possible      | <input type="checkbox"/> primary wound closure not possible                     |
|                       | <input type="checkbox"/> intact papillae                                | <input type="checkbox"/> impaired, missing papillae                             |
|                       | <input type="checkbox"/> adequate keratinised mucosa                    | <input type="checkbox"/> inadequate keratinised mucosa                          |
|                       | <input checked="" type="checkbox"/> uneventful                          |   |
| Implantation          | <input type="checkbox"/> simultaneously with bone augmentation (1 step) | <input checked="" type="checkbox"/> successively to bone augmentation (2 steps) |

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©Geistlich Pharma AG